

## JEDI Program/JDI01

### Compound Screening for SARS-CoV-2 Proteins Using MST/Dianthus

nsp12

nanoDSF / MST/ TRIC measurements

May 27, 2021



## ■ nanoDSF:

- nsp12 shows good thermal stability ( $T_m$  of  $\sim 47^\circ\text{C}$ ).
- Reducing agents and detergents have no significant impact on protein thermal stability, and the protein tolerates up to 3% DMSO.
- However, none of the tool cpds (GS-443902 trisodium = active Remdesivir metabolite), Ribavirin (ICN-1229) and Favipiravir (MCE)) nor ATP show significant effects on protein thermal stability.

## ■ TRIC (Dianthus):

- Nsp12 was successfully labelled with RED-Tris-NTA dye.
- However, none of the tool cpds (GS-443902 trisodium = active Remdesivir metabolite), Ribavirin (ICN-1229) and Favipiravir (MCE)) show binding to RED-Tris-NTA labelled nsp12.

## ■ Labelled MST:

- For comparison, experiments were also run on the NT.115 instrument and are in line with TRIC and nanoDSF results.
- Nsp12 alone is insufficient for binding of ATP or tool compounds. Thus, a TRIC/MST assay cannot be established with nsp12.

# nanoDSF


nsp12 (DVT1, PC13929-1)

# nanoDSF Assay Conditions

<b>Protein</b>	0.1 mg/ml nsp12 (DVT1, PC13929-1)
<b>Assay Buffer</b>	8 Buffers (see table)
<b>Instrument</b>	Prometheus NT.48
<b>Capillary type</b>	nanoDSF Standard Capillaries
<b>Measurement parameters</b>	LED Power: 50 % Temperature ramp: 2°C/min

	Conditions
Buffer 1	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub>
Buffer 2	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>0.005 % Tween20</b>
Buffer 3	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>0.05 % Tween20</b>
Buffer 4	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>0.05 % Pluronic</b>
Buffer 5	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>0.1 % Pluronic</b>
Buffer 6	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>1 mM TCEP</b>
Buffer 7	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>2.5 mM DTT</b>
Buffer 8	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , <b>2.5 mM GSH</b>

# Effect of reducing agents and detergents on nsp12 (DVT1, PC13929-1) thermal stability



Condition	$\bar{\theta} T_m$ [°C] <sup>1</sup>	s.d. [°C]	$\bar{\theta} T_{agg}$ [°C] <sup>1</sup>	s.d. [°C]	Analysis mode
Buffer 1	47.6	0.2	45.8	0.1	ratio
Buffer 2	47.1	0.3	44.1	0.1	ratio
Buffer 3	46.6	0.0	43.7	0.0	ratio
Buffer 4	47.2	0.1	45.4	0.1	ratio
Buffer 5	47.2	0.1	45.0	0.1	ratio
Buffer 6	47.3	0.2	45.6	0.1	ratio
Buffer 7	47.2	0.1	45.3	0.6	ratio
Buffer 8	47.3	0.1	45.9	0.2	ratio


<sup>1</sup> determined in duplicate

- No significant impact of detergents or reducing agents on protein thermal stability was observed

# nanoDSF Assay Conditions

<b>Protein</b>	0.1 mg/ml nsp12 (DVT1, PC13929-1)
<b>Assay Buffer</b>	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.005 % Tween20 DMSO: 0 - 5 %
<b>Instrument</b>	Prometheus NT.48
<b>Capillary type</b>	nanoDSF Standard Capillaries
<b>Measurement parameters</b>	LED Power: 50 % Temperature ramp: 2°C/min

# Impact of DMSO on nsp12 (DVT1, PC13929-1) thermal stability



— 0% DMSO  
— 1% DMSO  
— 2% DMSO  
— 3% DMSO  
— 4% DMSO  
— 5% DMSO

DMSO	T <sub>m</sub> [°C]	s.d. [°C]	ΔT <sub>m</sub> [°C] <sup>2</sup>	Ø T <sub>agg</sub> [°C] <sup>1</sup>	s.d. [°C]	Analysis mode
-	46.8	0.1	-	43.6	1.6	ratio
1%	46.6	0.2	- 0.2	43.6	0.3	ratio
2%	46.0	0.1	- 0.8	39.6	4.6	ratio
3%	46.3	0.1	- 0.5	43.7	0.2	ratio
4%	45.6	0.2	- 1.2	42.9	0.1	ratio
5%	45.5	0.2	- 1.3	42.2	0.4	ratio



- Increasing DMSO concentrations destabilize the protein.
- A concentration of 2.5 % was selected for further measurements.

# nanoDSF Assay Conditions

<b>Protein</b>	0.1 mg/ml nsp12 (DVT1, PC13929-1)
<b>Assay Buffer</b>	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.005 % Tween20 DMSO: 2.5 %
<b>compounds</b>	ATP at 100 $\mu$ M and 500 $\mu$ M GS-443902 trisodium (JDI-880) at 10 $\mu$ M and 100 $\mu$ M Ribavirin (ICN-1229) (JDI-879) at 10 $\mu$ M and 100 $\mu$ M Favipiravir (MCE) (JDI-881) at 10 $\mu$ M and 100 $\mu$ M
<b>Instrument</b>	Prometheus NT.48
<b>Capillary type</b>	nanoDSF Standard Capillaries
<b>Measurement parameters</b>	LED Power: 50 % Temperature ramp: 2°C/min



# Effect of compounds on nsp12 (DVT1, PC13929-1) thermal stability



Titrant	$\emptyset T_m$ [°C] <sup>1</sup>	s.d. [°C]	$\Delta T_m$ [°C] <sup>2</sup>	$\emptyset T_{agg}$ [°C]	s.d. [°C]	Analysis mode
Protein alone	46.9	0.1	-	44.7	0.2	ratio
2.5 % DMSO	46.4	0.1	- 0.5	44.3	0.5	ratio
500 $\mu$ M ATP	47.3	0.1	0.4	44.8	0.5	ratio
100 $\mu$ M ATP	47.1	0.2	0.2	45.3	0.3	ratio
100 $\mu$ M GS-443902 trisodium	46.5	0.2	- 0.4	44.3	0.1	ratio
10 $\mu$ M GS-443902 trisodium	46.5	0.2	- 0.4	44.5	0.0	ratio
100 $\mu$ M Ribavirin (ICN-1229)	45.8	0.1	- 1.1	43.1	0.0	ratio
10 $\mu$ M Ribavirin (ICN-1229)	46.6	0.4	- 0.3	44.7	0.1	ratio
100 $\mu$ M Favipiravir (MCE)	47.4	0.0	0.5	45.6	0.2	ratio
10 $\mu$ M Favipiravir (MCE)	47.1	0.1	0.2	44.7	0.5	ratio

<sup>1</sup> determined in duplicate

<sup>2</sup> referenced to

- Neither ATP nor any of the tool compounds has a significant impact on protein thermal stability.

# nanoDSF Summary: cpds effects on nsp12 (DVT1, PC13929-1)

Effect of 3 compounds and ATP on the thermal stability of nsp12 (DVT1, PC13929-1) analyzed by nanoDSF

Target name	Titrant / Condition	$\emptyset T_m$ [°C]	s.d. [°C]	$\Delta T_m$ [°C]	$\emptyset T_{agg}$ [°C]	s.d. [°C]	Analysis mode
nsp12	Protein alone	46.9	0.1	-	44.7	0.2	ratio
	2.5 % DMSO	46.4	0.1	- 0.5	44.3	0.5	ratio
	500 $\mu$ M ATP	47.3	0.1	0.4	44.8	0.5	ratio
	100 $\mu$ M ATP	47.1	0.2	0.2	45.3	0.3	ratio
	100 $\mu$ M GS-443902 trisodium	46.5	0.2	- 0.4	44.3	0.1	ratio
	10 $\mu$ M GS-443902 trisodium	46.5	0.2	- 0.4	44.5	0.0	ratio
	100 $\mu$ M Ribavirin (ICN-1229)	45.8	0.1	- 1.1	43.1	0.0	ratio
	10 $\mu$ M Ribavirin (ICN-1229)	46.6	0.4	- 0.3	44.7	0.1	ratio
	100 $\mu$ M Favipiravir (MCE)	47.4	0.0	0.5	45.6	0.2	ratio
	10 $\mu$ M Favipiravir (MCE)	47.1	0.1	0.2	44.7	0.5	ratio


# TRIC (Dianthus)

nsp12 (DVT1, PC13929-1)


# TRIC Assay Conditions

<b>Fluor. Molecule</b>	25 nM RED-tris-NTA 1 <sup>st</sup> gen. 25 nM RED-tris-NTA 2 <sup>nd</sup> gen.	
<b>Fluorophore</b>	RED-tris-NTA 1 <sup>st</sup> and 2 <sup>nd</sup> gen.	
<b>Instrument</b>	Dianthus NT.23PicoDuo	
<b>Measurement parameter</b>	LED Power: 11 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Singlicates	
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.05% Tween20 Buffer 2: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.5 mM ATP, 0.05% Tween20	
<b>Titrant</b>	nsp12 (DVT1, PC13929-1)	1000 – 0.49 nM (12 conc.)

# RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. vs. nsp12 (DVT1, PC13929-1) at Buffer 1



RED-tris-NTA 1<sup>st</sup> gen.




RED-tris-NTA 2<sup>nd</sup> gen.


Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	Lower confidence [M]	Upper confidence [M]	$\Delta F_{norm}$ [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	3.3E-09	1.8E-09	5.9E-09	46.6	24.2	5	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	RED-tris-NTA 2 <sup>nd</sup> gen.	nsp12	1.0E-10	1.1E-14	8.8E-07	55.9	14.4	5	Buffer 1

- RED-tris-NTA 1<sup>st</sup> gen. binds nsp12 with a determined  $K_D$  of 3.3 nM.
- RED-tris-NTA 2<sup>nd</sup> gen. binds nsp12 with a determined  $K_D$  of 0.1 nM.
- Both dyes can be used to non-covalently label nsp12.

# RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. vs. nsp12 (DVT1, PC13929-1) at Buffer 2



RED-tris-NTA 1<sup>st</sup> gen.



RED-tris-NTA 2<sup>nd</sup> gen.


Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	1.0E-10	-	-	40.5	43.6	5	Buffer 2
RED-tris-NTA 2 <sup>nd</sup> gen.	RED-tris-NTA 2 <sup>nd</sup> gen.	nsp12	1.7E-09	1.4E-09	2.0E-09	65.1	122.3	5	Buffer 2

- RED-tris-NTA 1<sup>st</sup> gen. binds nsp12 with a determined K<sub>D</sub> of 1.0 nM.
- RED-tris-NTA 2<sup>nd</sup> gen. binds nsp12 with a determined K<sub>D</sub> of 1.7 nM.
- Both dyes can be used to non-covalently label nsp12.



# TRIC Assay Conditions

<b>Fluor. Molecule</b>	25 nM nsp12 (DVT1, PC13929-1)		
<b>Fluorophore</b>	RED-tris-NTA 1 <sup>st</sup> and 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	25 nM protein / 12.5 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 23 % (nano detector) TRIC settings: 1 - 5 - 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.05% Tween20 Buffer 2: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.5 mM ATP, 0.05% Tween20		
<b>Titrant</b>	GS-443902 trisodium (active Remdesivir metabolite)	JDI-880	20 µM – 0.11 nM (12 conc.)


# RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. vs. nsp12 (DVT1, PC13929-1) in Buffer 1



RED-tris-NTA 1<sup>st</sup> gen.



RED-tris-NTA 2<sup>nd</sup> gen.




Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	-	5	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	-	5	Buffer 1


- RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. labelled nsp12 does not bind GS-443902 trisodium in buffer 1.



# RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. vs. nsp12 (DVT1, PC13929-1) in Buffer 2



RED-tris-NTA 1<sup>st</sup> gen.



RED-tris-NTA 2<sup>nd</sup> gen.


Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	-	5	Buffer 2
RED-tris-NTA 2 <sup>nd</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	-	5	Buffer 2

- RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. labelled nsp12 does not bind GS-443902 trisodium in buffer 2.

# TRIC Assay Conditions

<b>Fluor. Molecule</b>	25 nM nsp12 (DVT1, PC13929-1)		
<b>Fluorophore</b>	RED-tris-NTA 1 <sup>st</sup>		
<b>Labelling conditions</b>	25 nM protein / 12.5 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 20 % (nano detector) TRIC settings: 1 - 5 - 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.05% Tween20 DMSO: 2.5%		
<b>Titrant</b>	GS-443902 trisodium	JDI-880	200 µM – 1.13 nM (12 conc.)
	Ribavirin (ICN-1229)	JDI-879	
	Favipiravir (MCE)	JDI-881	


# RED-tris-NTA 1<sup>st</sup> labelled nsp12 (DVT1, PC13929-1) vs. GS-443902 trisodium



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	-	5	

- RED-tris-NTA labelled nsp12 does not bind GS-443902 trisodium.


# RED-tris-NTA 1<sup>st</sup> labelled nsp12 (DVT1, PC13929-1) vs. Ribavirin (ICN-1229)



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	Ribavirin (ICN-1229)	-	-	-	-	-	5	

- RED-tris-NTA labelled nsp12 does not bind Ribavirin (ICN-1229).

# RED-tris-NTA 1<sup>st</sup> labelled nsp12 (DVT1, PC13929-1) vs. Favipiravir (MCE)



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	Favipiravir (MCE)	-	-	-	-	-	5	

- RED-tris-NTA labelled nsp12 does not bind Favipiravir (MCE).

# Labelled MST


nsp12 (DVT1, PC13929-1)

# MST labelled assay conditions


<b>Fluor. Molecule</b>	25 nM nsp12 (DVT1, PC13929-1)		
<b>Fluorophore</b>	RED-tris-NTA 1 <sup>st</sup> and 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	25 nM protein / 12.5 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Monolith NT.115 (03)		
<b>Capillary type</b>	Monolith™ NT.115 Series MST Premium Coated Capillaries		
<b>Measurement parameter</b>	LED Power: 80 % MST Power: 40 % MST settings: 3 – 20 – 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicate		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.05% Tween20 Buffer 2: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.5 mM ATP, 0.05% Tween20		
<b>Titrant</b>	GS-443902 trisodium	JDI-880	20 µM – 0.11 nM (12 conc.)

# RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. labelled nsp12 (DVT1, PC13929-1) vs. GS-443902 trisodium in buffer 1

RED-tris-NTA 1st gen.



RED-tris-NTA 2nd gen.




Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	K <sub>D</sub> Confidence [M]	ΔFnorm [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	20	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	20	Buffer 1

- RED-tris-NTA 1<sup>st</sup> gen. labeled nsp12 does not bind GS-443902 trisodium.
- RED-tris-NTA 2<sup>nd</sup> gen. labeled nsp12 does not bind GS-443902 trisodium.




# RED-tris-NTA 1<sup>st</sup> and 2<sup>nd</sup> gen. labelled nsp12 (DVT1, PC13929-1) vs. GS-443902 trisodium in buffer 2

RED-tris-NTA 1st gen.



RED-tris-NTA 2nd gen.




Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	K <sub>D</sub> Confidence [M]	ΔFnorm [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	20	Buffer 2
RED-tris-NTA 2 <sup>nd</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	20	Buffer 2

- RED-tris-NTA 1<sup>st</sup> gen. labelled nsp12 does not bind GS-443902 trisodium.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled nsp12 does not bind GS-443902 trisodium.

# MST labelled assay conditions

<b>Fluor. Molecule</b>	25 nM nsp12 (DVT1, PC13929-1)		
<b>Fluorophore</b>	RED-tris-NTA 1 <sup>st</sup> and 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	25 nM protein / 12.5 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Monolith NT.115 (03)		
<b>Capillary type</b>	Monolith™ NT.115 Series MST Premium Coated Capillaries		
<b>Measurement parameter</b>	LED Power: 100 % MST Power: 40 % MST settings: 3 – 20 – 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicate		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl <sub>2</sub> , 2.5 mM DTT, 0.05% Tween20 DMSO: 2.5%		
<b>Titrant</b>	GS-443902 trisodium	JDI-880	200 µM – 1.13 nM (12 conc.)
	Ribavirin (ICN-1229)	JDI-879	
	Favipiravir (MCE)	JDI-881	


# [RED-tris-NTA 1<sup>st</sup> labelled nsp12 (DVT1, PC13929-1) vs. GS-443902 trisodium



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	K <sub>D</sub> Confidence [M]	ΔFnorm [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	GS-443902 trisodium	-	-	-	-	20	

- RED-tris-NTA labelled nsp12 does not bind GS-443902 trisodium.


# [RED-tris-NTA 1<sup>st</sup> labelled nsp12 (DVT1, PC13929-1) vs. Ribavirin (ICN-1229)



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	K <sub>D</sub> Confidence [M]	ΔFnorm [%o]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	Ribavirin (ICN-1229)	-	-	-	-	20	

- RED-tris-NTA labelled nsp12 does not bind Ribavirin (ICN-1229).

# RED-tris-NTA 1<sup>st</sup> labelled nsp12 (DVT1, PC13929-1) vs. Favipiravir (MCE)



Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	$K_D$ Confidence [M]	$\Delta F_{norm}$ [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 1 <sup>st</sup> gen.	nsp12	Favipiravir (MCE)	-	-	-	-	20	

- RED-tris-NTA labelled nsp12 does not bind Favipiravir (MCE).

## Next steps

- As neither nsp12 nor the RdRp complex allow biophysical characterization of compound interactions using TRIC (Dianthus), we recommend to stop the project

**Crelux GmbH** a WuXi AppTec company

Dr. Saskia Villinger

Senior Scientist and Deputy Head of Biophysics & Screening

Am Klopferspitz 19a

82152 Martinsried

Germany

[Saskia.Villinger@wuxiapptec.com](mailto:Saskia.Villinger@wuxiapptec.com)

[www.crelux.com](http://www.crelux.com)

[www.wuxiapptec.com](http://www.wuxiapptec.com)

