

JEDI Program/JDI01

Compound Screening for SARS-CoV-2 Proteins Using MST/Dianthus

Nsp12 - TRIC Based 12-Point K_D Validation

September 16, 2021



■ TRIC (Dianthus)

12-pt K_D validation on RED-TRIS NTA 2nd gen. labelled nsp12 was performed with the 8 most promising hits from 8-pt screen. Compounds were measured from 250 μ M, showing the following results:

- Two weak binders were identified: nsp12S5746081 and nsp1295277234
- Three compounds are non-binders
- Three compounds do not show binding, but lead to aggregation

TRIC (Dianthus)

Nsp12 (DVT1, PC13929-1)

vs. 8 Hits

TRIC Screening Conditions:

Fluor. Molecule	25 nM nsp12 (DVT1, PC13929-1)		
Fluorophore	12.5 nM RED-TRIS NTA 2 nd gen.		
Labelling conditions	Labelling buffer: 20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl ₂ , 2.5 mM DTT, 0.005% Tween 25 nM protein / 12.5 nM dye Incubation time: 60 min Centrifugation: 10 min at 15000g		
Preparation steps	Dilute Fragments with the Integra Viaflow Head (12-Point dilution series) Spotting of 500 nl compound with Mosquito into Dianthus Assay plate Buffer/labelled protein were added using Integra Pipettes and mixed with BioShake 3000 elm		
Instrument	Dianthus NT.23PicoDuo		
Measurement parameter	LED Power: 24 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicate		
Assay buffer	20 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl ₂ , 2.5 mM DTT, 0.005% Tween DMSO: 2.5 %		
Titrant	8 Hits	See list	250 µM – 1.41 nM (12 concentrations)
	Suramin		100 µM – 0.56 nM (12 concentrations)

Compounds

JDI-38

JDI-44

JDI-49

JDI-80


JDI-88

JDI-92

JDI-137

JDI-147


RED-TRIS NTA 2nd gen. labelled nsp12 vs. Suramin (Detector 1)



Fluorophore	Fluor. Molecule	Titrant	K _D [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	Suramin	6.3E-07	4.2E-07	9.4E-07	-28.2	16.4	5	

Suramin binds to RED-TRIS NTA 2nd gen. labelled nsp12 with estimated K_D value of 629 nM and stable ΔFnorm and signal-to-noise over the entire assay plate.



RED-TRIS NTA 2nd gen. labelled nsp12 vs. Suramin (Detector 2)



Fluorophore	Fluor. Molecule	Titrant	K _D [M]	Lower confidence [M]	Upper confidence [M]	ΔF_{norm} [‰]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	Suramin	5.3E-07	4.1E-07	6.9E-07	-26.7	24.8	5	

Suramin binds to RED-TRIS NTA 2nd gen. labelled nsp12 with estimated K_D value of 528 nM and stable ΔF_{norm} and signal-to-noise over the entire assay plate.



RED-TRIS NTA 2nd gen. labelled nsp12 vs. weak binder



Fluorophore	Fluor. Molecule	Titrant	K _D [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	nsp12S5746081 (JDI-49)	>3.9E-05	-10.9	9.2	5	Weak Binder, no saturation
RED-TRIS NTA 2nd gen.	Nsp12	nsp1295277234 (JDI-88)	>2.0E-04	-8.1	9.9	5	Weak Binder, no saturation

- RED-TRIS NTA 2nd gen. labelled nsp12 weakly binds nsp1242190400 with an estimated K_D > 40 μM without reaching saturation.
- RED-TRIS NTA 2nd gen. labelled nsp12 weakly binds nps5-45651217 with an estimated K_D > 200 μM without reaching saturation.


RED-TRIS NTA 2nd gen. labelled nsp12 vs. Aggregation



Fluorophore	Fluor. Molecule	Titrant	Category	K _D [M]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	N_mli9974624 (JDI-120)	Aggregation	-	-	5	Aggregation
RED-TRIS NTA 2nd gen.	Nsp12	nsp12S6174387 (JDI-132)	Aggregation	-	-	5	Aggregation

- RED-TRIS NTA 2nd gen. labelled nsp12 does not bind nsp1242190400, aggregation detected
- RED-TRIS NTA 2nd gen. labelled nsp12 does not bind S_mli69415177, aggregation detected



RED-TRIS NTA 2nd gen. labelled nsp12 vs. Aggregation



Fluorophore	Fluor. Molecule	Titrant	Category	K _D [M]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	nsp12S2245233 (JDI-137)	Aggregation	-	-	5	Aggregation

- RED-TRIS NTA 2nd gen. labelled nsp12 does not bind nsp12S2245233, aggregation detected


RED-TRIS NTA 2nd gen. labelled nsp12 vs. Non-Binder



Fluorophore	Fluor. Molecule	Titrant	K _D [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	nps5-45651217 (JDI-44)	-	-	-	5	ΔFnorm is insignificant
RED-TRIS NTA 2nd gen.	Nsp12	Nsp1205917961 (JDI-92)	-	-	-	5	ΔFnorm is insignificant

- RED-TRIS NTA 2nd gen. labelled nsp12 does not bind nps5-45651217, ΔFnorm is insignificant
- RED-TRIS NTA 2nd gen. labelled nsp12 does not bind Nsp1205917961, ΔFnorm is insignificant

RED-TRIS NTA 2nd gen. labelled nsp12 vs. Non-Binder



Fluorophore	Fluor. Molecule	Titrant	K _D [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	nsp1270326070 (JDI-147)	-	-	-	5	Insufficient saturation

- RED-TRIS NTA 2nd gen. labelled nsp12 does not bind nsp1270326070, insufficient saturation

TRIC summary

Fluorophore	Fluor. Molecule	Titrant	Category	K _D [M]	Signal / Noise	TRIC On [s]	Comment
RED-TRIS NTA 2nd gen.	Nsp12	nsp12S5746081 (JDI-49)	>3.9E-05	-10.9	9.2	5	Weak Binder, no saturation
RED-TRIS NTA 2nd gen.	Nsp12	nsp1295277234 (JDI-88)	>2.0E-04	-8.1	9.9	5	Weak Binder, no saturation
RED-TRIS NTA 2nd gen.	Nsp12	N_mli9974624 (JDI-120)	Aggregation	-	-	5	Aggregation
RED-TRIS NTA 2nd gen.	Nsp12	nsp12S6174387 (JDI-132)	Aggregation	-	-	5	Aggregation
RED-TRIS NTA 2nd gen.	Nsp12	nsp12S2245233 (JDI-137)	Aggregation	-	-	5	Aggregation
RED-TRIS NTA 2nd gen.	Nsp12	nps5-45651217 (JDI-44)	-	-	-	5	ΔFnorm is insignificant
RED-TRIS NTA 2nd gen.	Nsp12	Nsp1205917961 (JDI-92)	-	-	-	5	ΔFnorm is insignificant
RED-TRIS NTA 2nd gen.	Nsp12	nsp1270326070 (JDI-147)	-	-	-	5	Insufficient saturation

Next steps

- Discuss the data and next steps

Crelux GmbH a WuXi AppTec company

Dr. Saskia Villinger

Senior Scientist and Deputy Head of Biophysics & Screening

Am Klopferspitz 19a

82152 Martinsried

Germany

Saskia.Villinger@wuxiapptec.com

www.crelux.com

www.wuxiapptec.com

