

## JEDI Program/JDI01

### Compound Screening for SARS-CoV-2 Proteins Using MST/Dianthus

#### SARS-CoV-2 Spike protein

#### nanoDSF / MST/ TRIC measurements

May 28, 2021



## ■ nanoDSF:

- VHH 12 (nanobody against Spike) shows good thermal stability and a single melting transition.
- However, the nanobody was sensitive to buffer exchange: about 50% of protein was lost during buffer exchange.

## ■ TRIC:

- Spike was labelled according to the previously established RED-Tris-NTA labeling protocol.
- No binding of the nanobody was observed with or without buffer exchange (high protein loss during buffer exchange).
- Binding of ACE2 was observed with a determined  $K_D$  of 64 nM. However, some aggregation was observed as well.

## ■ Labelled MST:

- For comparison, binding of ACE2 was measured on the NT.115 instrument and observed with a determined  $K_D$  of 5 nM. This affinity is in line with literature and previous in-house experience.

➤ Spike TRIC/MST assay setup may be feasible with ACE2, but requires further assay optimization.


# nanoDSF

VHH 12 (EEH1, PD14993-1) nanobody for Spike

# nanoDSF Assay Conditions

<b>Protein</b>	0.1 mg/ml VHH 12 (EEH1, PD14993-1)
<b>Assay Buffer</b>	PBS
<b>Instrument</b>	Prometheus NT.48
<b>Capillary type</b>	nanoDSF Standard Capillaries
<b>Measurement parameters</b>	LED Power: 100 % Temperature ramp: 2°C/min

# Thermal stability of the Spike-directed nanobody



Condition	Ø T <sub>m</sub> [°C] <sup>1</sup>	s.d. [°C]	Analysis mode
0.1 mg/ml VHH 12	61.0	0.1	ratio


<sup>1</sup> determined in duplicate

- VHH 12 (nanobody against Spike) shows good thermal stability and a single melting transition.

# nanoDSF Assay Conditions

<b>Protein</b>	0.1 mg/ml VHH 12 (EEH1, PD14993-1)
<b>Assay Buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.05 % Tween20 Buffer 2: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.05 % Pluronic
<b>Instrument</b>	Prometheus NT.48
<b>Capillary type</b>	nanoDSF Standard Capillaries
<b>Measurement parameters</b>	LED Power: 100 % Temperature ramp: 2°C/min

# Effect of buffer exchange with different buffer on VHH 12 (EEH1, PD14993-1) thermal stability



Condition	Ø T <sub>m</sub> [°C] <sup>1</sup>	s.d. [°C]	Analysis mode
Buffer 1 w/o buffer exchange	63.0	0.1	ratio
Buffer 1 with buffer exchange	62.8	0.0	ratio
Buffer 2 w/o buffer exchange	60.2	0.1	ratio
Buffer 2 with buffer exchange	60.5	0.0	ratio

<sup>1</sup> determined in duplicate

- During buffer exchange about 50% of nanobody was lost
- Buffer exchange has no significant impact on the protein thermal stability of the remaining protein after buffer exchange
- The nanobody tolerates Tween + PEG-8000 slightly better than Pluronic.

# TRIC (Dianthus)



SARS-CoV-2 Spike (DYG4)



# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM S (SARS-CoV-2) (DYG4, PD14787-1)		
<b>Fluorophore</b>	25 nM RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 50 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 20 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.05% Tween Buffer 2: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.05% Pluronic Buffer 3: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% Pluronic Buffer 4: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.1% Pluronic		
<b>Titrant</b>	VHH 12 (nanobody) With buffer exchange (~50% protein loss)	EEH1 (PD14993-1)	20 MM – 0.11 nM, 12 conc., 1:3



# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike vs. VHH 12



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	DYG4	VHH 12	-	-	-	-	-	5	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	DYG4	VHH 12	-	-	-	-	-	5	Buffer 2

- RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike does not bind VHH 12 in buffer 1.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike does not bind VHH 12 in buffer 2.

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike vs. VHH 12





Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	DYG4	VHH 12	-	-	-	-	-	5	Buffer 3
RED-tris-NTA 2 <sup>nd</sup> gen.	DYG4	VHH 12	1.2E-06	2.7E-07	5.4E-5	-3.6	5.0	5	Buffer 4, low amplitude and S/N

- RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike does not bind VHH 12 in buffer 3.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike binds VHH 12 with a determined K<sub>D</sub> of 1.2 μM in buffer 4. However, ΔFnorm and signal-to-noise ratio are low.

# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM S (SARS-CoV-2) (DYG4, PD14787-1)		
<b>Fluorophore</b>	25 nM RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 50 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 18 % (nano detector) TRIC settings: 1 - 5 - 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: PBS, 0.05% Tween Buffer 2: PBS, 0.1% PEG 8000, 0.05% Tween		
<b>Titrant</b>	VHH 12 (nanobody) w/o buffer exchange (storage buffer = PBS)	EEH1 (PD14993-1)	20 µM – 0.11 nM, 12 conc., 1:3

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike vs. VHH 12




Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	DYG4	VHH 12	-	-	-	-	-	5	Buffer 1, without exchange
RED-tris-NTA 2 <sup>nd</sup> gen.	DYG4	VHH 12	-	-	-	-	-	5	Buffer 2, without exchange

- RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike does not bind VHH 12 in buffer 1.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike does not bind VHH 12 in buffer 2.

# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM S (SARS-CoV-2) (DYG4, PD14787-1)		
<b>Fluorophore</b>	25 nM RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 50 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 16% (nano detector) TRIC settings: 1 - 5 - 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Pluronic, 0.1% PEG-8000		
<b>Titrant</b>	ACE2	DYF3 (PD14701-1)	500 nM – 0.122 nM, 12 conc., 1:2

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike (DYG4) vs. ACE2



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-Tris-NTA (2. gen.)	Spike	ACE2	6.4E-08	-	-	12.1	9.6	1.5	aggregation

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Spike binds ACE2 with a determined K<sub>D</sub> of 64 nM. However, some aggregation was observed.

# Labelled MST


SARS-CoV-2 Spike (DYG4)



# MST labelled assay conditions

<b>Fluor. Molecule</b>	50 nM S (SARS-CoV-2) (DYG4, PD14787-1)		
<b>Fluorophore</b>	25 nM RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 50 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Monolith NT.115 (03)		
<b>Capillary type</b>	Monolith™ NT.115 Series MST Premium Coated Capillaries		
<b>Measurement parameter</b>	LED Power: 80 % MST Power: 40 % MST settings: 3 – 15 – 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicate		
<b>Assay buffer</b>	20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Pluronic, 0.1% PEG-8000		
<b>Titrant</b>	ACE2	DYF3 (PD14701-1)	500 nM – 0.122 nM, 12 conc., 1:2

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Spike vs. ACE2



Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	$K_D$ Confidence [M]	$\Delta F_{\text{norm}}$ [%]	Signal / Noise	MST on [s]	Comment
RED-Tris-NTA (2. gen.)	Spike	ACE2	4.8E-09	1.3E-09 - 1.76E-08	20	17.6	2.5	

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Spike binds ACE2 with a determined  $K_D$  of 4.8 nM on the NT.115 instrument, in line with literature and previous in-house experience.

# Next steps

- Discuss nanobody QC
- Further TRIC/MST assay optimization (additional Spike and tag-less ACE2 ordered, awaiting protein receipt)

**Crelux GmbH** a WuXi AppTec company

Dr. Saskia Villinger

Senior Scientist and Deputy Head of Biophysics & Screening

Am Klopferspitz 19a

82152 Martinsried

Germany

[Saskia.Villinger@wuxiapptec.com](mailto:Saskia.Villinger@wuxiapptec.com)

[www.crelux.com](http://www.crelux.com)

[www.wuxiapptec.com](http://www.wuxiapptec.com)

