

JEDI

JDI01

SARS-CoV-2 Nucleocapsid protein

MST/ TRIC measurements

August 5, 2021



- **TRIC (Dianthus) and Labelled MST (NT.115):**
  - Binding of **PolyT hexamers** to **Nucleocapsid** was tested on the Dianthus and Monolith NT.115 in buffers with different salt concentrations.
    - At 150 mM NaCl, **weak binding without saturation** was observed on both instruments. **Aggregation** was observed on the Dianthus.
    - At lower salt concentrations (10 and 50 mM), no binding and aggregation was observed on both instruments.
  - **Nanobodies** were prepared with dialysis (instead of gel filtration), and binding of nanobodies to **Nucleocapsid** was tested and optimized on Dianthus/NT.115.
    - On the Dianthus, **binding** of both nanobodies to Nucleocapsid was observed with high affinity. However, strong **aggregation** was observed and  $K_D$  values could not be fitted.
    - Buffer optimization was attempted in PBS buffers (w/o buffer exchange) and Hepes buffers (w/ dialysis). No improvement was observed, **aggregation** could not be prevented.
    - On the NT.115 instrument (MST), better data quality with less aggregation was observed:
      - RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid binds **VHH H3-3** with a determined  $K_D$  of 10 nM.
      - RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid binds **VHH E4-3** with a determined  $K_D$  of 93 nM.
  - **We recommend to use MST with one nanobody as tool compound for the compound screening.**



# TRIC (Dianthus)

SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)

# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 37 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, <b>150 mM NaCl</b> , 0.1% PEG 8000, 0.05% Tween20 Buffer 2: 20 mM HEPES pH 7.5, <b>50 mM NaCl</b> , 0.1% PEG 8000, 0.05% Tween20 Buffer 3: 10 mM Na-phosphate pH 7.5, <b>10 mM NaCl</b> , 0.1% PEG 8000, 0.05% Tween20 DMSO: 0%		
<b>Titrant</b>	PolyT	JDI-884	45 µM – 22 nM (12 conc.), 1:2


# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. PolyT



Fluorophore	Fluor. Molecule	Titant	$K_D$ [M]	Lower confidence [M]	Upper confidence [M]	$\Delta F_{\text{norm}}$ [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	PolyT	>1.9E-05	-	-	19.2	16.1	5	Buffer 1, aggregation, no saturation
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	PolyT	-	-	-	-	-	5	Buffer 2, aggregation

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid binds PolyT with an estimated  $K_D > 19 \mu\text{M}$  in buffer 1 without reaching saturation. In addition, aggregation was observed.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind PolyT in buffer 2, aggregation was observed.

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. PolyT




Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	PolyT	-	-	-	-	-	5	Buffer 3, aggregation

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind PolyT in buffer 3, strong aggregation was observed.


# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 28 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.05% Tween20 DMSO: 0%		
<b>Titrant</b>	VHH H3-3 (nanobody against CTD) VHH E4-3 (NTD nanobody) Dialyzed into: 20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Tween, 0.1% PEG 8000	EEG1 (PD14991-1) (stored at 4°C) EEF1 (PD14989-1) (stored at 4°C)	1 µM – 0.49 nM (12 conc.)

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3 and VHH E4-3



VHH H3-3  
(EEG1)



VHH E4-3  
(EEF1)

Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	aggregation
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH E4-3	-	-	-	-	-	5	aggregation



- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid shows binding to VHH H3-3 and VHH E4-3. However, due to aggregation, no K<sub>D</sub> value can be determined.



# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 41 % (nano detector) TRIC settings: 1 - 5 - 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: <b>PBS</b> , 0.1% PEG 8000, 0.05% Tween20 Buffer 2: <b>PBS</b> , 5% Glycerol, 0.05% Tween20 Buffer 3: <b>PBS</b> , 0.1% PEG 8000, 5% Glycerol, 0.05% Tween20 Buffer 4: <b>PBS</b> , 0.1% BSA, 0.05% Tween20 Buffer 5: <b>PBS</b> , 0.1% PEG 8000, 2mM DTT, 0.05% Tween20 DMSO: 0%		
<b>Titrant</b>	VHH H3-3 (nanobody against CTD)	EEG1 (PD14991-1) (stored at 4°C)	1 µM – 0.49 nM (12 conc.)



# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 2

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 1, aggregation was observed.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 2, aggregation was observed.


# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3



Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	Lower confidence [M]	Upper confidence [M]	$\Delta F_{\text{norm}}$ [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 3
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 4

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 3, aggregation was observed.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 4, aggregation was observed.

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3





Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 5

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 5, aggregation was observed.

# TRIC Assay Conditions

<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Dianthus NT.23PicoDuo		
<b>Measurement parameter</b>	LED Power: 28 % (nano detector) TRIC settings: 1 - 5 - 1 (s)      (initial fluorescence – MST on time – back-diffusion) Duplicates		
<b>Assay buffer</b>	Buffer 1: 20 mM <b>HEPES</b> pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.005% Tween20 Buffer 2: 20 mM <b>HEPES</b> pH 7.5, 150 mM NaCl, 5% Glycerol, 0.005% Tween20 Buffer 3: 20 mM <b>HEPES</b> pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 5% Glycerol, 0.005% Tween20 Buffer 4: 20 mM <b>HEPES</b> pH 7.5, 150 mM NaCl, 0.1% BSA, 0.005% Tween20 Buffer 5: 20 mM <b>HEPES</b> pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 2mM DTT, 0.005% Tween20 DMSO: 0%		
<b>Titrant</b>	VHH H3-3 (nanobody against CTD) Dialyzed into: 20 mM Hepes pH 7.5, 150 mM NaCl, 0.005% Tween	EEG1 (PD14991-1) (stored at 4°C)	1 µM – 0.49 nM (12 conc.)



# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 2

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 1, aggregation was observed.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 2, aggregation was observed.


# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	$\Delta F_{\text{norm}}$ [‰]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 3
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 4

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 3, aggregation was observed.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 4, aggregation was observed.

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	Lower confidence [M]	Upper confidence [M]	ΔFnorm [%]	Signal / Noise	TRIC On [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	-	5	Buffer 5

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 5, aggregation was observed.



# Labelled MST


SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)

# MST labelled assay conditions


<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Monolith NT.115 (03)		
<b>Capillary type</b>	Monolith™ NT.115 Series MST Premium Coated Capillaries		
<b>Measurement parameter</b>	LED Power: 100 % MST Power: 40 % MST settings: 3 – 20 – 1 (s)      (initial fluorescence – MST on time – back-diffusion) Singlicate and Duplicate		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.05% Tween20 Buffer 2: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.05% Tween20 Buffer 3: 10 mM Na-phosphate pH 7.5, 10 mM NaCl, 0.1% PEG 8000, 0.05% Tween20 DMSO: 0%		
<b>Titrant</b>	PolyT	JDI-884	45 µM – 22 nM (12 conc.), 1:2

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. PolyT

Buffer 1




Buffer 2



Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	$K_D$ Confidence [M]	$\Delta F_{\text{norm}}$ [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	PolyT	$> 2.4\text{E-}05$	-	17.3	11.3	10	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	PolyT	-	-	-	-	10	Buffer 2, aggregation

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid binds PolyT with an estimated  $K_D > 24 \mu\text{M}$  without reaching saturation in buffer 1.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind PolyT in buffer 2 and aggregation was observed.

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. PolyT



Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	$K_D$ Confidence [M]	$\Delta F_{norm}$ [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	PolyT	-	-	-	-	10	Buffer 3


- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind PolyT in buffer 3 and aggregation was observed.

# MST labelled assay conditions


<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Monolith NT.115 (03)		
<b>Capillary type</b>	Monolith™ NT.115 Series MST Premium Coated Capillaries		
<b>Measurement parameter</b>	LED Power: 80 % MST Power: 40 % MST settings: 3 – 20 – 1 (s) (initial fluorescence – MST on time – back-diffusion) Singlicate		
<b>Assay buffer</b>	20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.05% Tween20 DMSO: 0%		
<b>Titrant</b>	VHH H3-3 (nanobody against CTD) VHH E4-3 (NTD nanobody) Dialyzed into: 20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Tween, 0.1% PEG 8000	EEG1 (PD14991-1) (stored at 4°C) EEF1 (PD14989-1) (stored at 4°C)	1 µM – 0.49 nM (12 conc.)

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3 and VHH E4-3

VHH H3-3



VHH E4-3



Fluorophore	Fluor. Molecule	Titrant	$K_D$ [M]	$K_D$ Confidence [M]	$\Delta F_{norm}$ [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	1.0E-08	1.7E-09 – 6.2E-08	18.7	10.1	10	VHH H3-3
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH E4-3	9.3E-08	6.0E-08 – 1.5E-07	26.2	35.1	10	VHH E4-3


- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid binds VHH H3-3 with a determined  $K_D$  of 10 nM.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid binds VHH E4-3 with a determined  $K_D$  of 93 nM.
- Low amount of aggregation was observed for both nanobodies.

# MST labelled assay conditions


<b>Fluor. Molecule</b>	50 nM SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1)		
<b>Fluorophore</b>	RED-tris-NTA 2 <sup>nd</sup> gen.		
<b>Labelling conditions</b>	100 nM protein / 25 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g		
<b>Instrument</b>	Monolith NT.115 (03)		
<b>Capillary type</b>	Monolith™ NT.115 Series MST Premium Coated Capillaries		
<b>Measurement parameter</b>	LED Power: 80 % MST Power: 40 % MST settings: 3 – 20 – 1 (s)      (initial fluorescence – MST on time – back-diffusion) Singlicate and Duplicate		
<b>Assay buffer</b>	Buffer 1: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.05% Tween20 Buffer 3: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 5% Glycerol, 0.05% Tween20 Buffer 4: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% BSA, 0.05% Tween20 Buffer 5: 20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 2mM DTT, 0.05% Tween20 DMSO: 0%		
<b>Titrant</b>	VHH H3-3 (nanobody against CTD) Dialyzed into: 20 mM Hepes pH 7.5, 150 mM NaCl, 0.005% Tween	EEG1 (PD14991-1) (stored at 4°C)	1 µM – 0.49 nM (12 conc.)

# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3

Buffer 1



Buffer 3




Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	K <sub>D</sub> Confidence [M]	ΔFnorm [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	10	Buffer 1
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	10	Buffer 3

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 1.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 3.




# RED-tris-NTA 2<sup>nd</sup> gen. labelled SARS-CoV-2 Nucleocapsid protein (ECJ1, PD14785-1) vs. VHH H3-3

Buffer 4



Buffer 5



Fluorophore	Fluor. Molecule	Titrant	K <sub>D</sub> [M]	K <sub>D</sub> Confidence [M]	ΔFnorm [%]	Signal / Noise	MST on [s]	Comment
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	10	Buffer 4
RED-tris-NTA 2 <sup>nd</sup> gen.	Nucleocapsid	VHH H3-3	-	-	-	-	10	Buffer 5

- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 4.
- RED-tris-NTA 2<sup>nd</sup> gen. labelled Nucleocapsid does not bind VHH H3-3 in buffer 5.

# Next steps

- We recommend to test the **robustness of MST** in the best assay buffer (20 mM HEPES pH 7.5, 150 mM NaCl, 0.1% PEG 8000, 0.05% Tween20) using VHH H3-3 or VHH E4-3 nanobody as tool protein (prepared by dialysis).
  - If suitable, screen the compounds for Nucleocapsid binding using MST instead of Dianthus.
  - If necessary, perform additional MST buffer optimization.
- Potentially: Further Dianthus assay optimization, which we estimate to be very challenging with this protein batch

**Crelux GmbH** a WuXi AppTec company

Dr. Saskia Villinger

Senior Scientist and Deputy Head of Biophysics & Screening

Am Klopferspitz 19a

82152 Martinsried

Germany

[Saskia.Villinger@wuxiapptec.com](mailto:Saskia.Villinger@wuxiapptec.com)

[www.crelux.com](http://www.crelux.com)

[www.wuxiapptec.com](http://www.wuxiapptec.com)

