

JEDI Program/JDI01

Compound Screening for SARS-CoV-2 Proteins Using MST/Dianthus

SARS-CoV-2 Spike protein

TRIC (Dianthus), MST

July 23, 2021



■ TRIC/Dianthus:

- A pretest was performed to verify sufficient binding of RED-Tris-NTA dye to the Spike protein provided by Pasteur (DYG6).
 - Binding was observed with a K_D of 4 nM, thus allowing Spike labelling via the His-tag.
- The binding assay was tested using the Pasteur protein and two ACE2 batches on the Dianthus.
 - Binding of Fc-tagged ACE2 was observed with $K_D < 1$ nM. However, aggregation was also observed.
 - No binding of untagged ACE2 was observed.

■ Labelled MST:

- In parallel, the established assay was performed with the Pasteur Spike protein on the Monolith NT.115 instrument.
 - Binding of Fc-tagged ACE2 was observed with low nM affinity, in line with the K_D observed for the commercial Spike protein batch (DYG4).
 - Binding of untagged ACE2 was observed with estimated $K_D > 190$ nM without reaching saturation. Aggregation was also observed.
- Overall, labelled MST (NT.115 instrument) is feasible with Spike using Fc-tagged ACE2 as positive control.
- TRIC (Dianthus) assay optimization is very challenging and would require more efforts and protein.


TRIC (Dianthus)

SARS-CoV-2 Spike provided by Pasteur (DYG6)

TRIC Assay Conditions: Pretest

| | | | |
|------------------------------|--|-------------------|--------------------------|
| Fluor. Molecule | 12.5 nM RED-tris-NTA 2 nd gen. | | |
| Fluorophore | RED-tris-NTA 2 nd gen. | | |
| Instrument | Dianthus NT.23PicoDuo | | |
| Measurement parameter | LED Power: 23 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicates | | |
| Assay buffer | 20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Tween, 0.1% PEG-8k | | |
| Titrant | Pasteur Spike (SARS-CoV-2) | (DYG6, PD15372-1) | 350 – 0.17 nM (12 conc.) |

RED-tris-NTA 2nd vs. SARS-CoV-2 Spike provided by Pasteur (DYG6)




| Fluorophore | Fluor. Molecule | Titrant | EC ₅₀ [M] | Lower confidence [M] | Upper confidence [M] | Hill coeff. | ΔFnorm [‰] | Signal / Noise | TRIC On [s] | Comment |
|-----------------------------------|----------------------------------|---------------|----------------------|----------------------|----------------------|-------------|------------|----------------|-------------|-----------|
| RED-tris-NTA 2 nd gen. | RED-tris-NTA 2 nd gen | Pasteur Spike | 4.1E-09 | - | - | 2.6 | 14.0 | 14.0 | 1.5 | 2nd event |

- RED-tris-NTA 2nd gen. dye binds with low nM affinity to Pasteur Spike (EC₅₀ of 4 nM, K_D fit does not represent the data). At higher concentrations of Spike, a second event is observed that is difficult to interpret.


TRIC Assay Conditions

| | | | |
|------------------------------|---|------------------|--------------------------|
| Fluor. Molecule | 50 nM Pasteur Spike (SARS-CoV-2) (DYG6, PD15372-1) | | |
| Fluorophore | 25 nM RED-tris-NTA 2 nd gen. | | |
| Labelling conditions | 100 nM protein / 50 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g | | |
| Instrument | Dianthus NT.23PicoDuo | | |
| Measurement parameter | LED Power: 29 % (nano detector) TRIC settings: 1 - 5 - 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicates | | |
| Assay buffer | 20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Tween, 0.1% PEG-8k | | |
| Titrant | ACE2 Fc-tagged | DYF1 (PD13357-1) | 400 – 0.20 nM (12 conc.) |
| | ACE2 untagged | DYF4 (PD15147-1) | 500 – 0.24 nM (12 conc.) |
| | Preparation with dialysis tubes | | |

RED-tris-NTA 2nd gen. labelled SARS-CoV-2 Pasteur Spike (DYG6) vs. ACE2 (Fc-tagged, DYF1 and untagged, DYF4)



ACE2 Fc-tagged



ACE2 untagged

| Fluorophore | Fluor. Molecule | Titrant | K _D [M] | Lower confidence [M] | Upper confidence [M] | ΔFnorm [‰] | Signal / Noise | TRIC On [s] | Comment |
|-----------------------------------|----------------------|----------------|--------------------|----------------------|----------------------|------------|----------------|-------------|--------------------------------|
| RED-tris-NTA 2 nd gen. | Pasteur Spike (DYG6) | ACE2 Fc-tagged | < 1.0E-9 | - | - | 19.0 | 14.2 | 1.5 | strong binder, aggregation |
| RED-tris-NTA 2 nd gen. | Pasteur Spike (DYG6) | ACE2 untagged | - | - | - | - | - | 1.5 | No binding, strong aggregation |

- RED-tris-NTA 2nd gen. labelled Spike binds ACE2 Fc-tagged with a determined K_D < 1 nM. However, aggregation was observed at later TRIC on times.
- RED-tris-NTA 2nd gen. labelled Spike does not bind ACE2 untagged. Aggregation was observed as well.



Labelled MST

SARS-CoV-2 Spike provided by Pasteur (DYG6)

MST labelled assay conditions

| | | | |
|------------------------------|---|------------------|--------------------------|
| Fluor. Molecule | 50 nM Pasteur Spike (SARS-CoV-2) (DYG6, PD15372-1) | | |
| Fluorophore | 25 nM RED-tris-NTA 2 nd gen. | | |
| Labelling conditions | 100 nM protein / 50 nM dye Incubation time: 30 min Centrifugation: 10 min at 15000g | | |
| Instrument | Monolith NT.115 (03) | | |
| Capillary type | Monolith™ NT.115 Series MST Premium Coated Capillaries | | |
| Measurement parameter | LED Power: 50 % MST Power: 40 % MST settings: 3 – 15 – 1 (s) (initial fluorescence – MST on time – back-diffusion) Duplicate | | |
| Assay buffer | 20 mM Hepes pH 7.5, 150 mM NaCl, 0.05% Tween, 0.1% PEG-8k | | |
| Titrant | ACE2 Fc-tagged | DYF1 (PD13357-1) | 400 – 0.20 nM (12 conc.) |
| | ACE2 untagged | DYF4 (PD15147-1) | 500 – 0.24 nM (12 conc.) |
| | Preparation with dialysis tubes | | |



RED-tris-NTA 2nd gen. labelled SARS-CoV-2 Spike vs. ACE2 (Fc-tagged, DYF1)



| Fluorophore | Fluor. Molecule | Titrant | K_D [M] | K_D Confidence [M] | ΔF_{norm} [%] or ΔF | Signal / Noise | MST on [s] | Comment |
|-----------------------------------|----------------------|----------------|-----------|----------------------|-------------------------------------|----------------|----------------------|---------|
| RED-tris-NTA 2 nd gen. | Pasteur Spike (DYG6) | ACE2 Fc-tagged | 1.2E-09 | 4.0E-10 – 3.6E-09 | 28.6 | 37.2 | 1.5 | |
| RED-tris-NTA 2 nd gen. | Pasteur Spike (DYG6) | ACE2 untagged | 4.0E-09 | 3.0E-09 – 5.3E-09 | 81.1 | 19.5 | Initial Fluorescence | |

- RED-tris-NTA 2nd gen. labelled Spike binds ACE2 Fc-tagged with a determined K_D of 1.2 and 4 nM, using MST and initial fluorescence, respectively. Due to the high protein concentration in the assay setup, the K_D may be error-prone.

RED-tris-NTA 2nd gen. labelled SARS-CoV-2 Spike vs. ACE2 (untagged, DYF4)



| Fluorophore | Fluor. Molecule | Titrant | K_D [M] | K_D Confidence [M] | ΔF_{norm} [‰] or ΔF | Signal / Noise | MST on [s] | Comment |
|-----------------------------------|----------------------|---------------|-------------------|----------------------|--|----------------|----------------------|----------------------------|
| RED-tris-NTA 2 nd gen. | Pasteur Spike (DYG6) | ACE2 untagged | $>1.9\text{E-}07$ | - | - | 11.2 | 1.5 | No saturation, aggregation |
| RED-tris-NTA 2 nd gen. | Pasteur Spike (DYG6) | ACE2 untagged | - | - | - | - | Initial Fluorescence | ΔF insignificant |

- RED-tris-NTA 2nd gen. labelled Spike binds ACE2 untagged with estimated $K_D > 190 \text{ nM}$ without reaching saturation. Aggregation is observed.

Next steps

- Assay setup is ready with commercial Spike and Fc-tagged ACE2 (positive control) using labelled MST (NT.115):
 - Measure the compounds in 8 pt screening (singlicate)
- Potentially: further assay optimization using the Dianthus, which we estimate to be very challenging for this protein target

Crelux GmbH a WuXi AppTec company

Dr. Saskia Villinger

Senior Scientist and Deputy Head of Biophysics & Screening

Am Klopferspitz 19a

82152 Martinsried

Germany

Saskia.Villinger@wuxiapptec.com

www.crelux.com

www.wuxiapptec.com

